

HI·bebe

BT-200 Fetal Doppler

Pocket Size Doppler

for
**Pregnant Women
Midwives
Obstetricians**



HI·bebe is the first way you can meet your baby

- Compact and Light
- High-Sensitivity Doppler Probe
- Ergonomic Design
- Large-Display & Clean Sound
- Easy to Use
- Low Power Consumption
- Long Time to Continuous Use (6 hours)

HI·bebe Fetal Doppler



Specification

Display	LCD
	Low Battery Display
	Signal Quality Display
Heart Rate Range	50~240
Accuracy	±2%
Ultrasound Frequency	2MHz
Ultrasound Intensity	<10mW/cm ²
Speaker Output Power	1.2W (Built-in Loudspeaker)
Audio output for earphone	
Auto Shut-off	5 minutes
Auto Shut-off on signa	1 minute
Power	1.5V Battery X 2 (Type:AA)
Sensitivity	10~12 Weeks Onward
PC Interface	Sound Card
Battery Life	360 minutes
Accessories	Carrying Case
	Battery(1.5V AA × 2)
	Gel(60ml)



Sound Version Available!

- No LCD Display
- Other functions are same