




Antepartum Fetal Monitor Model BT-300



BT-300 Antepartum fetal monitor sets new standards for office environment monitoring. Convenient, portable and Easy to use, the Model BT-300 is the best antepartum fetal monitor that meets user' satisfaction with all of these distinctive features:

- Compact and light
- High-Sensitivity from the 1 MHz probe
- Simultaneous external twin FHR with dual pulsed doppler
- Ergonomic designed Event marker
- Economic with outstanding Two Years Warranty

Fetal Monitor Model BT-300



Specifications

FHR

Dual Pulsed Doppler
Ultrasound Frequency : 1 MHz
Intensity : 10 mW/cm² or less
FHR Range : 50 ~ 240 bpm
FHR Accuracy : $\pm 2\%$ of Range

FM

Dual Fetal Movement

UC

External Type
Frequency Response : DC ~ 0.5 Hz
Reference (zero) Control
Measurement Range : 0 ~ 99 units

Printer

Thermal Array Type
Print Speed : 1, 2, 3 cm/min
Paper Feeding Function

Display

7-Segment LED
3 Channels (FHR I, FHR II, UC)

Sound

Doppler Sound
Alarm sound
Information Sound

Set-up

Alarm Upper / Lower Limit Value
Print Speed, Print Contrast

Function

Mark Function
FHR II Offset Function
Auto Print Function

Power

Input : AC 100 ~ 240 V (50/60 Hz)
Output : DC 16 V, 2.8 A
Consumption : 80 VA

Dimension

191mm(H) X 187mm(L) X 201mm(D)

External Link

RS-232C

Standard Accessory

Ultrasound Doppler Probe 2ea
UC Probe 1ea
Event Marker Jack 1ea
Print Paper 2ea
Power Adaptor & Cord 1ea
Ultrasound Gel 1ea
Probe Belt 3ea
Operation Manual 1ea

Warranty

2 Years (probe excluded)

Distributors address

201, 2F, 239-15, Gasan-Dong, Geumcheon-Gu, Seoul, Korea
TEL 82-2-869-0642, 0643 FAX 82-2-862-0644
E-Mail ssrobber@hotmail.com

 **bistos**
Bistos Co., Ltd