

Fetal Monitor BT-350

BT-350 is amazing color LCD fetal monitor to check your fetal health

- 7" wide TFT color LCD display
- Up to 150 patients data saving
- Fetal heart sound play and record in PC
- Screen is rotating to 180°
- Multiple language support
- Quick guide display
- Various installation method



BT-350 LED



BT-350 LCD

Fetal Monitor BT-350

Specifications

FHR	• Dual Pulsed Doppler • Intensity : 10 mW / cm² or less • Accuracy : $\pm 2\%$ of Range	• Ultrasound Frequency : 0.985 MHz • FHR Range : 30 ~ 240 bpm FHR
FM	• Auto-detection of Dual Fetal Movement	
UC	• External Type • Reference (zero) Control	• Frequency Response : DC~0.5 HZ • Measurement Range : 0 ~ 99 units
Printer	• Thermal Array Type • Auto print : Off, 10, 20, 30, 40, 50, 60 min	• Print Speed : 1,2,3 cm/min, high speed • Paper Feeding Function
Display	• 3 Channels (FHR I, FHR II, UC) • 7" Wide TFT-Color LCD : 800(W) × 480(H) BT-350 LCD • Large 7-segment & volume level LED : BT-350 LED	
Sound	• Doppler Sound • Information Sound	• Alarm Sound
Set-up	• Alarm Upper / Lower Value • Print Speed	• Time & Date
Function	• Mark function • Auto Print function • Multi-language & quick guide function : BT-350 LCD	• FHR II Offset function • Fetal movement print function
Trend	• Data Saving for 450 hours(3hr / person) : BT-350 LCD	
Power	• Input : AC 100 ~ 240 V (50 / 60 Hz) • Battery Operation : continuously 2 hours	• Output : DC 18 V, 2.5 A Consumption
Dimension	• 80 mm(H) X 330 mm(L) X 280 mm(D)	
External Link	• Wire : RS-232 C	• Wireless : Bluetooth
Accessory Standard	• Ultrasound Doppler Probe 2 ea • Event Marker Jack 1ea • Power Adaptor & Cord 1 ea • Probe Belt 3 ea • AST(Acoustic Stimulator) 1 ea • Bluetooth Antenna	• UC Probe 1 ea • Print Paper 2 ea • Ultrasound Gel 1 ea • Operation Manual 1 ea • Battery
Optional		
Warranty	• 2 Years (Accessory excluded)	

PC Interface (central monitoring system)

BCM-330

- Remote control up to 8 monitors
- Fetal heart sound Pla and Record
- Unlimited Multi language support
- Wire (RS-232C) and Wireless (Bluetooth) connection



Notebook type



Stand type



Wall-mount type