

LEISEGANG® LM-900 Cryosurgery System

Efficient and Safe Freeze for
Cryotherapy Procedures



LEISEGANG® LM-900 Cryosurgery System, the benchmark in cervical care, offers an effective, efficient and safe freeze for cryotherapy procedures

Features and Benefits:

- Excellent thermal conductivity, designed for continuous usage
- Calibrated, easy-to-read pressure gauge ensures accurate freeze temperatures
- Insulated probe barrel, high-stress gas lines and solid brass manifold offer added safety
- Instant off and automatic gas line cleaning with remote on/off switch
- Solid, lightweight body design, comfort-shaped handle and finger-trigger control provide ease of operation with less fatigue
- Space-saving, compact console with integrated tip holder offers added convenience
- Wide choice of handcrafted, double-bonded, coin-silver cryotip probes for every lesion and location (cold soak and EO sterilization only)

Order as a system package and save:

- LM-900 with custom case
- Empty 20 lb gas cylinder (specify N₂O or CO₂) and mobile cart
- O-rings (10)
- Choice of four probes with cleaning plugs

Recommended OB/GYN probe tips:

- 50105 General Purpose
- 50108 General Purpose
- 50109 Exocervical
- 50110 Exo/Endocervical

LM-900 Probes




Microderm #50101



Soll Lash Epilating
#50102



Condyloma #50103



Derm Bevel #50104



General Purpose*
#50105



Large Derm Bevel
#50106



Curved Anorectal
#50107



General Purpose*
#50108



Exocervical* #50109



Exo/Endocervical*
#50110



Special Exo/Endo
#50111



Exo/Endocervical
#50112



Maxiderm #50113



Large Exocervical
#50114



Jumbo Exo/Endo
#50115

**Best Sellers*

To find out more or to purchase the LM-900 Cryosurgery System
call 800.243.2974.

CooperSurgical

CooperSurgical is the leading company dedicated to providing medical devices and procedure solutions that improve health care delivery to women regardless of clinical setting. Our company is fostering that position through expansion of its core businesses and introduction of advanced technology-based products which aid clinicians in the management and treatment of commonly seen conditions.

203.601.5200 | 800.243.2974 | coopersurgical.com | 95 Corporate Dr. | Trumbull, CT 06611